

## Epredia Single Use Disposal Cytofunnels

### Contact us for more pricing details!

#### The original, trusted, easy-to-use funnel

Designed for ease of use and efficiency, the single-use Epredia™ Shandon™ Cytofunnels™ are an ideal solution for today's busy laboratories.

A cost-effective and time-effective way to make quality thin-layer cell preparations. Designed for use with the Epredia™ Shandon™ Cytospin 4 Cytocentrifuge.

#### Cytofunnel general features

- Fit snugly against a specially marked Cytoslide
- Held in place with a stainless steel Cytoclip
- Cytoclips easily unlock to dispose of the sample chamber and to remove the slide
- Twelve assemblies fit into the Cytospin head
- Include caps and easy-to-use operational illustrations

#### EZ Cytofunnel general features

- For single use with all Epredia Cytospin cytocentrifuge applications
- Optimises cell yield
- Produces cost-effective, thin-layer preparations
- Reduces the risk of exposure to pathogenic samples
- Disposable to eliminate time-consuming decontamination and cleaning
- Decreases the possibility of cross-contamination
- Eliminates the need for metal Shandon™ Cytoclips™ slide clip

Description	Code
<b>Cytofunnel</b>	
Single Cytofunnel White with Cap 50pcs	EPKM5991040
Single Cytofunnel White with Cap 500pcs	EPKM1102548
Double Cytofunnel White with Cap 25pcs	EPKM5991039
Double Cytofunnel White with Cap 500pcs	EPKM1102547
<b>EZ Cytofunnel</b>	
EZ Cytofunnel White with Cap 40pcs	EPKMA78710003
EZ Cytofunnel White with Cap 200pcs	EPKMA78710020



Normal Cytofunnel



EZ CytoFunnel

POWERED BY SHANDON  MICROM  MENZEL-GLÄSER  Richard-Allan Scientific 

Contact us for more information  
and to place your orders

#### Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: sales.au@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 September 2021

Ref. Code:  
EPKM-CYTO

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.