

Tomy SX-700E Autoclave

Affordable Sterilising Solution.

This top loading autoclave incorporates a 79 litre stainless steel chamber with two baskets. With integrated (but removable) condensate bottle to trap condensate from chamber during venting.

- Runs on 240V 50Hz 15amp supply
- Includes a condensate bottle - no plumbing required!
- Mobile - can be easily moved
- 1 free preventative maintenance visit (metro area's only)
- Digital temperature and time display with analogue gauge for pressure
- Simple lid opening using either hand or foot
- Cooling fan to reduce run time
- Five sterilisation profiles to select from (in addition to normal sterilisation)
- Three operating temperature ranges - warming 45 to 95°C, Heating 45 to 104°C and Sterilising 105 to 135°C
- Safety devices include - water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection and safety valve



**Save Over
\$7,000!**

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Autoclave 79L with Two Baskets	TOMYSX-700E	\$17,086	\$9,975

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



VALID UNTIL:
30 September 2021

Ref. Code:
TOMY-SX700E

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.