

BioEcho 50 Reaction Kits

Introductory Offer!

Next Generation kits DNA/RNA extraction kits, providing more convenience, quality and speed.

By using the principle of reverse chromatography, BioEcho extraction methods ensure DNA remains untouched and any unwanted sample-derived contaminant is bound and eliminated by the column. Potential inhibitory process reagents are not in use and thus do not have to be removed.

The result: pure DNA by only one centrifugation step.

After just one single spin, the purified high molecular weight RNA/DNA is eluted, while all debris and cellular compartments, inhibitors, small fragments (< 50 bp), nucleotides, and salts are completely removed.

30% Off



By moving to a BioEcho extraction kit you are supporting:

- environmentally friendly science: no harsh extraction chemicals
- sustainability: with 70% less plastics than standard Silica kits

Code	Description	List	Sale
BIOE010-011-050	EchoLUTION Blood DNA HiYield Kit	\$260	\$182
BIOE010-001-050	EchoLUTION Blood DNA Micro Kit	\$240	\$168
BIOE010-006-050	EchoLUTION Cell Culture DNA Kit	\$260	\$182
BIOE010-003-050	EchoLUTION Plant DNA Kit	\$300	\$210
BIOE010-002-050	EchoLUTION Tissue DNA Micro Kit	\$310	\$217
BIOE012-001-050	EchoLUTION Viral RNA/DNA Swab Kit	\$330	\$231
BIOE012-002-050	EchoLUTION Viral RNA/DNA Swab Kit Plus	\$370	\$259
BIOE020-002-040-050	EchoCLEAN Organic Solvent DNA CleanUp Kit	\$240	\$168
BIOE020-002-050-050	EchoCLEAN RNA CleanUp Kit	\$240	\$168



Contact us for more information and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2022

Ref. Code:
BIOE-50s

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.