

## Epredia Cytospin™ 4 Cyto-centrifuge

### The ultimate, affordable system for thin-layer cell preparation.

**Purchase a Cytospin 4 Cyto-centrifuge and receive 1 x Clipped Cytofunnel Starter Kit free of charge (valued at \$1600).**

The Epredia Cytospin™ 4 has been synonymous with cyto-centrifugation for over 40 years.

As 'Gold Standard' in cyto-centrifugation, more than 16,000 laboratories worldwide use this simple, reliable instrument. It provides economical thin-layer preparations from any liquid matrix, especially hypocellular fluids such as spinal fluids and urines.

The Cytospin 4 utilises carefully controlled centrifugation to separate and deposit a thin layer of cells on slides while maintaining cellular integrity.

This proven technology delivers excellent representation of all cell types present in homogenous liquid samples.

**Free kit includes: 50 Single White Cytofunnels, 25 Double Cytofunnels, Cytoslides, 12 Cytoclips.**



**Free  
Starter Kit!**

Description	Code
Cytospin 4 Cyto-centrifuge	EPRNA78300002
FREE Clipped Cytofunnel Starter Kit	EPRN3120110

POWERED BY

SHANDON

MICROM

MENZEL-GLÄSER

Richard-Allan  
Scientific

Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 September 2021

Ref. Code:  
EPRN-KIT

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.