

## IKA KS4000i Control Incubator with AS 4000.1 Universal Attachment

**Only 2 available at this special price!**

**Hurry stocks won't last.**

Innovative incubator shaker design allowing unattended operation in a temperature-controlled environment.

- Large LED display for speed and time settings
- Controls with antimicrobial coating for reduction of bacteria
- Integrated PID temperature control (use of two PT 1000 temperature sensors)
- Junction box in the workspace for connection of an additional temperature sensor e.g. PT 1000.60
- Electronic temperature and speed control
- Electronic timer switch: 0 - 999 h (set by the minute or by the hour)
- Optional built-in cooler for connection to an external cooling unit e.g. RC 2 basic / control
- Unit switches off automatically if disturbed
- Unit stops automatically when hood is lifted
- Collecting tray with drain hose on rear of unit
- Simple operation
- All functions can be controlled and documented with labworldsoft®

**Special  
Offer!**



Description	Code	Price (excl. GST)	Sale Price (excl. GST)
IKA KS 4000 i Control Special Offer – Includes AS 4000.1 Universal Attachment	IKAA3510001	\$13,285	<b>\$9,995</b>

### **IKA RCT Basic on sale for only \$850!**

**The must have for all laboratories!**

Contact Bio-Strategy today.

Description	Code
IKA RCT Basic Hotplate Stirrer	IKAA20002620

The improvement of the bestseller:  
Now with new technology for more capacity.

New stronger motor for a higher range of speeds.

New additional temperature control mode for faster heating of medium.



Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 September 2021  
or until sold out

Ref. Code:  
IKAA-KS4

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.