

## Labconco SteamScrubber and FlaskScrubber

### Contact Bio-Strategy for more information.

The all new GlassScrubber for general labware and FlaskScrubber for narrow-necked glassware, Labconco laboratory glassware washers have specialised features to meet the laboratory's demand for superior cleaning and convenience.

- Underbench or freestanding models
- Up to 6 kw heating
- Save 40% water usage when compared to hand-washing
- CleanWorks™ OS for intelligent wash programming

SteamScrubber Glassware Washers are designed primarily to hold general labware such as beakers, utensils, test tubes, Petri dishes and BOD bottles. Upper and lower standard open racks (included) accommodate a variety of optional inserts.

FlaskScrubber Glassware Washers are designed for narrow-neck glassware such as volumetric and Erlenmeyer flasks. A lower spindle rack with 36 spindles (included) uses direct injection to clean labware inside and out. Optional specialised inserts to wash pipettes, BOD bottles, petri dishes, culture tubes, utensils, etc. are available.



**NEW!**

Contact us for more information  
and to place your orders

#### **Bio-Strategy Pty Ltd.**

T: 1800 008 453  
E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
[shop.bio-strategy.com](http://shop.bio-strategy.com)

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Rev. July 2021

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.