

Miele PG8583 CD Washer/Disinfector

Save \$1,879 on this package deal!

Free installation and commissioning (city metro areas only within 100km of GPO), free A102 upper basket and 2 AK12 inserts, free detergent and neutraliser (total savings: \$1,879).

This underbench /standalone model combines a powerful pump with integrated heating elements for safe washing cycles up to 93°C. The DryPlus active HEPA filtered drying ensures loads are dry and ready to use.

- 16 programs including new "Mini" program for lightly soiled glassware
- Conductivity monitoring – proof that load is "clean"
- Configured for washing wide neck and narrow neck glassware
- Dispensing pumps for liquid detergent and liquid neutraliser
- Compact: 900mm x 700mm x 835mm W x D x H (820mm H with lid removed)
- Pullout drawer for housing chemical canisters (detergent and neutraliser)
- Stainless steel interior/exterior for durability and cleanliness
- ADP model including pump for AD water and waterproofing kit
- Many other options and accessories to select from
- Spray arm and spray pressure monitoring

Package Deal!



Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Miele PG8583 CD Package 1 Includes: A102 Upper Basket with 2 x AK12 Inserts, A150 Lower Module Basket with 2 x A301/5 Modules (3 x 6 = 18 Nozzles/Module), and Waterproofing Kit	MIELPG8583CD	\$21,825	\$19,946

Contact us for more information and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



VALID UNTIL:
30 September 2021

Ref. Code:
MIEL-PG8583CD

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.