

SOURCE

For more information and product codes please go to www.bio-strategy.com/promotions

Save!

analytikjena
An Endress+Hauser Company

Analytik Jena SpeedMill PLUS

Modern preparation of samples with the Analytik Jena SpeedMill PLUS

Achieving high quality results starts with careful sample preparation including efficient and reproducible sample homogenization.

- Innovative double action technology
- Passive sample cooling to -80°C
- Compact design / minimal footprint
- 2-year warranty



Description	Code	Price	Sale
SpeedMill PLUS 230 V	ANJE845-00007-2	\$14,795	\$11,390

[View here for more details](#)

Introductory Offer!

GESIM

GeSim BioScaffold Printer BS 3.3 Prime

Bioprinting made easy

Tool head with three Z-axes

- Tools included: pneumatic cartridges for RT and high/low temperature nozzle set, 50ml syringe extruder, UV pen (many more tools can be added)
- Substrate tray
- Size: 740x570mm² (WxD, including space for cables)
- Tray size: 350x250mm²



[View here for more details](#)

Usually \$130,000

Introductory Offer:
\$100,000

In Stock!

TOMY

Tomy SX-700E Toploading Autoclave

Supplied with two stainless steel baskets and a 2-year warranty!

The Tomy SX-700E top loading autoclave is powered from a 240V 50Hz 15amp supply and does not require plumbing so can be moved to where it is required!

- Integrated (removable) condensate bottle to trap condensate from chamber during venting
- Simple lid opening using hand or foot
- Cooling fan to reduce run time
- 1 **Free** Preventative Maintenance visit (metro areas only)



Description	Code	Price	Sale
Autoclave, toploading, 79 litre	TOMYSX-700E	\$15,514	\$9,995

[View here for more details](#)

Free Accessory!

PHCbi

PHCbi Pharmaceutical Refrigerator

For a limited time, FREE OF CHARGE - Battery housing and battery accessory for power failure alarm

Add Quality of Cold to your cold chain with PHCbi MPR models.

- Superior temperature uniformity assures safety of your products
- Temperature independent auto defrost (no devastating temperature spikes!)
- Proven reliability you can depend on!

In Stock!



Description	Code
PHCbi 1364 L Double Door Pharmaceutical Refrigerator	SANYMPR-1412-PE
PHCbi 684 L Single Door Pharmaceutical Refrigerator	SANYMPR-722-PE
Battery Mounting Box for Power Failure Alarm (MPR-722-PE or MPR-1412-PE models) - FOC	SANYMPR-48B1-PW

[View here for more information](#)

SOURCE

For more information and product codes please go to www.bio-strategy.com/promotions

Special Offer!



ILC Dover SoloLAB

Introductory offer: \$3,000

The soloLAB is a pop-up, benchtop isolation unit which is super quick to assemble (it really does just pop-up) that offers containment for many applications. Utility sleeves allow gas lines to be added to control the environment within.



The soloLAB provides excellent value for money to any labs looking to contain a variety of processes or substances.

[View here for more details](#)

Free Printer!



Priorclave QCS 150L Top-Loading Electrically Heated Autoclave

Free Printer (valued at \$2,000)

The unit is suitable for a variety of loads including fermentation / reaction vessels. Supplied with three loading baskets.

- Media Warming Facility
- Single handed lever action door mechanism
- 10-year door and chamber warranty
- Elements in chamber for ease of servicing



Description	Code	Price	Sale
Priorclave QCS 150L Top-Loading Electrically Heated Autoclave	PRIOPS/QCS/EV150	\$32,392	\$28,795

[View here for more details](#)

100% Guarantee!



Bethyl ELISA Kits

High-performing ELISAs you can count on

Every Bethyl product ships with a 100% guarantee; choose Bethyl with confidence.

Bethyl Laboratories are proud to offer ELISA kits for the detection of human IgG and IgM antibodies against SARS-CoV-2 in human serum or plasma.



Description	Code
Human IgG ELISA Kit	BETHE88-104

[View here for more details](#)

10% Off!



CYTIVA Vivaspin™ Protein Concentrator Spin Columns

Top quality products on sale, order today!

Vivaspin™ ultrafiltration spin columns are designed for fast, non-denaturing concentration and desalting of proteins and other biological samples.

- Minimal protein adsorption – thanks to the Vivaspin polycarbonate tube design



Description	Code	Price	Sale
Concentrator MWCO 3000 VivaSpin 500 [25 item]	GEHE28-9322-18	\$132	\$119

[View here for more details](#)

20% Off!



Reichert Handheld Refractometers

Versatile, pocket-sized digital devices

From the versatile, pocket-sized digital devices for performance and value, to the durable and rugged analog, these instruments put powerful and accurate performance in the palm of your hand.



Description	Code	Price	Sale
RHINO™ BRIX50 REFRACTOMETER	REIC137531L0	\$450	\$360
BRIX/RI-CHEK	REIC13940000	\$680	\$544
AR200 W/ IR COMM. PKG.	REIC13952000	\$3,485	\$2,788

[View here for more details](#)

15% Off!



IKA - Discounted across our entire Product Range

IKA - Designed to work perfectly

Proven technology for mixing, heating, distilling and crushing applications. Order now and save!

Items Include:

- T 25 Digital ULTRA-TURRAX package
- Vortex 3
- Eurostar 60
- MultiDrive Control
- A11 Basic



(not available in conjunction with other special offers)

[View here for more details](#)

Package Deal!

Miele**Miele PG8583 Washer/Disinfector**

Free installation and commissioning*, free detergent and neutraliser (total savings: \$742)

Package also includes A102 upper basket with 2 x AK12 inserts, A150 lower module basket with 2 x A301/5 modules (3 x 6 = 18 nozzles/module).

Reduce the risk of injury through manual washing and improve safety in your wash area!



*(city metro areas only within 100km of GPO)

Description	Code	Price	Sale
Miele PG8583 Washer/Disinfector Complete Package 1	MIELPG-8583	\$13,802	\$12,400

[View here for more details](#)

Price Drop!

epredia**Epredia™ Histoplast PE Paraffin Wax****Resealable Bags, 1kg**

Designed to suit all histological applications, including accommodations for rapid infiltration, extra-thin sectioning, and hard/fibrous tissue

- Made of paraffin and plastic polymers
- Melting Point: 56 - 57°C
- Good for general infiltration and embedding procedures



Description	Code	Price	Sale
Histoplast PE Paraffin Wax (8 x 1kg)	EPKM8330	\$92.35	\$45

[View here for more details](#)

Stock Clearance!

eppendorf**Eppendorf Combitips Advanced® Clearance****Smart, Safe, and Efficient**

Large numbers of users worldwide value the Multipipette / Combitip system, an efficient tool for quick and precise filling of plates or large series of tubes as used in nucleic acid purification kits.

- 9 volume sizes from 0.1ml to 50ml
- 4 purity grades including Eppendorf Sterile, Quality, Biopur® and Forensic DNA Grade



Description	Code	Price	Sale
COMBITIP ADVANCED 0,1ML WHITE [100 items]	EPPE0030089.405	\$165	\$155

[View here for more details](#)

Contact Us Today!

CORNING**Corning® CoolCell®****Timesaving, Cost Effective, Easy Handling**

Corning CoolCell provides alcohol-free freezing at a rate of -1°C/min, ideal for cell cryopreservation.

- Rapidly returns to room temperature when removed from the freezer
- Lower cost of use than alcohol-based devices
- Ease of use
- Simple, consistent way to standardise controlled-rate freezing



Unless otherwise specified, all products are for research use only.

[View here for more details](#)

Save!

Techno Plas
SERVICE | INNOVATION | RESULTS**TechnoPlas 70ml Universal Flat Bottom Container (Screw Cap)**

Specimen containers for the collection of biological specimens.

They are also ideal for storage and water sampling.

- 3 start thread for operator ease of cap removal
- Leak proof qualities



Description	Code	Price	Sale
70ml Flat Bottom Container Unlabelled Yellow Cap [550 items]	TCHPP5744UUUV	\$87.20	\$75
70ml Flat Bottom Container Labelled Yellow Cap [550 items]	TCHPP5744ULUV	\$91.40	\$80

[View here for more details](#)

Introductory Offer!

PALL Laboratory**Pall Minimate EVO Tangential Flow Filtration System**

Streamlines laboratory-scale concentration, desalting, and buffer exchange processes.

System components are designed for easy assembly and disassembly. Tubing and fitting replacement is simple and quick. Slip nuts for luer locks eliminate tubing kinks.



Description	Code	Price	Sale
Minimate Evo Tangential Flow Filtration System	PALLOAPMPUNV	\$8,768	\$8,100

[View here for more details](#)

SOURCE

For more information and product codes please go to www.bio-strategy.com/promotions

Trade In!



Update Your Lab for 2022: Reader Trade-In Deals

Receive up to 35% discount*

Trade in any Microplate Reader to a new Molecular Devices Reader
Tells us your requirements and we can create a trade in quote tailor made just for you.

- Any make or model
- Working or not



Save up to
35% on a
new reader

** Trade In Offer Excludes - SpectraMax L, SpectraMax ABS and SpectraMax ABS Plus readers*

[View here for more details](#)

Demo Clearance!



Demo Stock Clearance

Molecular Devices IXM-Confocal System

A unique confocal imaging solution capable of imaging more than a million wells a week.

- Transmitted Light
- 3 disk geometries
- Various objectives, and filters

Available now!

Demo Stock
Save Over
60%



Contact us for more information

[View here for more details](#)

Great Discount!



UV/VIS Plate Reader for \$19,990

Special - not to be repeated!

Purchase the ABS PLUS reader for only \$19,990 (usually \$32,000)

- 96-well microplate, 384-well microplates and cuvettes
- 190-1000 nm
- Endpoint and kinetic and spectrum read modes
- 4 research licenses for SoftMax Pro 7.1




[View here for more details](#)

Enter to Win a Food Hamper delivered before the Holidays

Complete the below brainteaser for a chance to win a Food Hamper to share with your work mates.

"A scientist puts bacteria in a petri dish at exactly noon. Every minute the bacteria divides into two and doubles in size. At exactly 1pm the petri dish is full. At what time was the dish half full?"



[Register your answer here](#)

Entries Close 30th November 2021

Answer and Winner will be announced in our December eNewsletter

Contact us for more information and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)



[Buy Online](#)



[View Promotions](#)

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.