

## Eppendorf Combitips Advanced® Clearance

### Smart, Safe, and Efficient

#### An efficient tool for quick and precise filling of plates or large series of tubes as used in nucleic acid purification kits.

Thanks to their revolutionary design, Combitips advanced dispenser tips provide you with many benefits for work in a modern laboratory. With an integrated piston, they function according to the positive displacement principle and ensure that the correct volume of liquid is dispensed regardless of the physical properties of the liquid you work with (e.g., viscosity, volatility, density, temperature). You can therefore achieve high-precision dispensing, independent of the flow characteristics of the liquid, and safely work with radioactive or toxic materials due to the hermetically sealed piston that prevents aerosol contamination.

#### Key features:

- Available in 9 volume sizes from 0.1ml to 50ml
- Available in 4 purity grades including Eppendorf Sterile, Quality, Biopur® and Forensic DNA Grade
- The elongated tips of the Combitips advanced make it possible to empty all standard lab tubes entirely
- Individual color-coding makes it easy to quickly identify the desired Combitips advanced

| Description                     | Code            | List            | Sale         |
|---------------------------------|-----------------|-----------------|--------------|
| <b>COMBITIP ADVANCED</b>        |                 |                 |              |
| 0,1ML WHITE [100 ITEMS]         | EPPE0030089.405 |                 |              |
| 0,5ML PURPLE [100 ITEMS]        | EPPE0030089.421 |                 |              |
| 1ML YELLOW [100 ITEMS]          | EPPE0030089.430 |                 |              |
| 2,5ML GREEN [100 ITEMS]         | EPPE0030089.448 | <b>\$165.47</b> | <b>\$155</b> |
| 5ML BLUE [100 ITEMS]            | EPPE0030089.456 |                 |              |
| 10ML ORANGE [100 ITEMS]         | EPPE0030089.464 |                 |              |
| 25ML RED [100 ITEMS]            | EPPE0030089.472 | <b>\$195.31</b> | <b>\$185</b> |
| 50ML L-GREY [100 ITEMS]         | EPPE0030089.480 | <b>\$215.20</b> | <b>\$205</b> |
| 0,1ML WHITE BIOPUR [100 ITEMS]  | EPPE0030089.618 |                 |              |
| 0,2ML L-BLUE BIOPUR [100 ITEMS] | EPPE0030089.626 |                 |              |
| 0,5ML PURPLE BIOPUR [100 ITEMS] | EPPE0030089.634 |                 |              |
| 1,0ML YELLOW BIOPUR [100 ITEMS] | EPPE0030089.642 | <b>\$236.55</b> | <b>\$220</b> |
| 2,5ML GREEN BIOPUR [100 ITEMS]  | EPPE0030089.650 |                 |              |
| 5ML BLUE BIOPUR [100 ITEMS]     | EPPE0030089.669 |                 |              |
| 10ML ORANGE BIOPUR [100 ITEMS]  | EPPE0030089.677 |                 |              |
| 25ML RED BIOPUR [100 ITEMS]     | EPPE0030089.685 | <b>\$277.39</b> | <b>\$220</b> |
| 50ML L-GREY BIOPUR [100 ITEMS]  | EPPE0030089.693 | <b>\$298.51</b> | <b>\$280</b> |

**SAVE**



Contact us for more information  
and to place your orders

#### Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: sales.au@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
31 December 2021

Ref. Code:  
EPPE-COMBITIPS

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.