



## Molecular Devices ABS Plus Reader

### UV/VIS Plate Reader - One off Opportunity

Purchase the ABS PLUS reader for  
only \$19,990 (usually \$32,000)

The SpectraMax ABS Plus absorbance microplate readers provide the flexibility, sensitivity, and convenience for a wide range of assays such as ELISAs, microbial growth, and protein quantitation.

This reader features an 8-channel read head for fast reads and a subset feature the patented PathCheck™ Sensor.

The temperature independent PathCheck sensor automatically normalizes and ensures accurate absorbance values by detecting pipetting errors.

- Free installation, training and application support
- 1-year warranty
- 96-well microplate, 384-well microplates and cuvettes
- 190-1000 nm
- Endpoint and kinetic and spectrum read modes
- 4 research licenses for SoftMax Pro 7.1

**Bonus  
Offer!**



Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453

E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

[shop.bio-strategy.com](http://shop.bio-strategy.com)