

## Tomy SX-700E Autoclave

### Top-loading 79 litre including extras

Supplied with two stainless steel baskets and a 2-year warranty!

The Tomy SX-700E top-loading autoclave is powered from a 240V 50Hz 15amp supply and does not require plumbing so can be moved to where it is required!

- Integrated (removable) condensate bottle to trap condensate from chamber during venting.
- Simple lid opening using either hand or foot
- Cooling fan to reduce run time
- 1 **Free** Preventative Maintenance visit  
*(city metro area's only within 100km of GPO)*
- Three operating temperature ranges:  
Warming 45 to 95°C, Heating 45 to 104°C and Sterilizing 105 to 135°C
- Safety devices include:
  - water level sensor
  - current leakage breaker
  - lid interlock
  - overheat prevention
  - over pressure prevention
  - open temperature sensor detection & safety valve
- Digital temperature & time display with analogue gauge for pressure.
- Five sterilization profiles to select from (in addition to normal sterilization)

**In Stock!**



Description	Code	List	Sale
Autoclave, toploading, 79 litre	TOMYSX-700E	\$15,514	<b>\$9,995</b>

Contact us for more information and to place your orders

#### Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: sales.au@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
31 December 2021

Ref. Code:  
TOMY-SX700E

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.