

Bio Echo - New Supplier

Bio-Strategy now supplying BioEcho across Australia

BioEcho has a completely new approach for DNA/RNA extraction and analysis.

Next Generation kits for DNA/RNA extraction and analysis. Based on the principle of reverse chromatography, the DNA remains untouched, any unwanted sample-derived contaminant is bound and eliminated by the column, potential inhibitory process reagents are not in use and, thus, don't have to be removed.

The result: pure DNA by only one centrifugation step, getting an inhibitor-free PCR template.

The complete BioEcho portfolio consists of:

- EchoLUTION DNA extraction kits
- EchoCLEAN DNA cleanup kits
- EchoSAFE buffers and reagents
- ReDUCE Sample concentration
- Lab tools



BioEcho's products are suited for multiple areas of biotechnology, including genomics, transcriptomics, next generation sequencing, biobanking, forensics, pathology, veterinary, food research, as well as industrial, medical and academic research.

Contact us or [visit our website](#) today for more information

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com