

Corning Squat Beakers

PYREX® Beaker with Spout

For 100 years, Corning has developed special glass for use in both chemical and life science laboratories, including PYREX® glass. Made from Type 1, Class A low expansion borosilicate glass, PYREX glassware has become the accepted standard in chemistry labs across the globe.

Low Form, Double Scale, Graduated

Manufactured with:

- Uniform wall thickness,
- Optimum balance between thermal shock resistance and mechanical strength
- Double graduated metric scale to indicate approximate content
- Extra large marking spot

| Description | Code | Price | Sale |
|-----------------------------------|----------------|---------------|------------|
| BEAKER SQUAT GLASS 50ML [1 ITEM] | CORN1000-50_U | \$4.20 | \$2 |
| BEAKER SQUAT GLASS 250ML [1 ITEM] | CORN1000-250_U | \$5.80 | \$3 |
| BEAKER SQUAT GLASS 600ML [1 ITEM] | CORN1000-600_U | \$7.70 | \$4 |

Heavy Duty Beakers

Specifically designed to offer:

- Best mechanical strength under harsh conditions, such as mechanized washing operations
- The 250 through 4000 ml beakers have a double graduated metric scale to indicate their approximate content.
- All sizes have an extra large marking spot

| Description | Code | Price | Sale |
|--|----------------|----------------|------------|
| BEAKER SQUAT HEAVY DUTY GLASS 150ml [1 ITEM] | CORN1003-150_U | \$4.10 | \$3 |
| BEAKER GRADUATED HEAVY DUTY - 250ml [1 ITEM] | CORN1003-250 | \$7.00 | \$3 |
| BEAKER SQUAT HEAVY DUTY GLASS 400ml [1 ITEM] | CORN1003-400_U | \$9.10 | \$4 |
| BEAKER SQUAT HEAVY DUTY GLASS 600ml [1 ITEM] | CORN1003-600_U | \$10.70 | \$5 |
| BEAKER SQUAT HEAVY DUTY GLASS 1L [1 ITEM] | CORN1003-1L_U | \$20.40 | \$9 |

**Stock
Clearance!**



Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2021

Ref. Code:
CORN_Beaker

only while stock lasts

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.