

## Miele PG-8593 Washer / Disinfector

### Washer Disinfector with DryPlus HEPA Filtered Air Drying

**Free external DOS K85 flex dispensing/dosing pump for detergent or neutraliser plus 5L container of liquid detergent and 5L container of neutraliser**

Laboratory Underbench/Free-Standing Washer with stainless steel interior combining powerful pump with integrated heating elements enabling safe washing cycles up to 93°C .

- 15 programs including new "Mini" program for lightly soiled glassware
- Internal dispensing pump for liquid detergent or liquid neutraliser
- Requires optional upper and lower baskets with inserts for a full working system
- Spray arm and spray pressure monitoring with alerts to advise you if a blockage or pressure deviation has occurred
- Compact: 600mm x 600mm x 850mm (820mm with lid removed)
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- Stainless steel interior/exterior for durability and cleanliness
- AD model for AD water and waterproofing kit
- Many other options and accessories to select from

**Package Deal**  
**Save \$959!**



Description	Code	List	Sale
Miele PG-8583 Washer / Disinfector Package 1	MIELPG-8593 PK	\$16,341	\$15,382

Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453  
E: sales.au@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com