

Tomy SX-700E Toploading Autoclave

Take advantage of 2021 pricing while stocks last!

Supplied with two stainless steel baskets and a 2-year warranty!

The Tomy SX-700E top loading autoclave is powered from a 240V 50Hz 15amp supply and does not require plumbing so can be moved to where it is required!.

- Integrated (removable) condensate bottle to trap condensate from chamber during venting
- Simple lid opening using either hand or foot
- Cooling fan to reduce run time
- Three operating temperature ranges: **Warming** 45 to 95°C, **Heating** 45 to 104°C and **Sterilizing** 105 to 135°C
- Safety devices include: water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection & safety valve
- Digital temperature & time display with analogue gauge for pressure
- Five sterilization profiles to select from (in addition to normal sterilization)
- 1 **Free** Preventative Maintenance visit (*metro area's only*)

**In
Stock!**



Description	Code	Price	Sale
Mobile 79 Litre Toploading Autoclave	TOMYSX-700E	\$15,514	\$9,995

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2022

Ref. Code:
TOMY-SX-700E

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.