

## Cytiva Whatman

### Grade 1 Qualitative Filter Papers

The most widely used filter paper for routine laboratory applications, qualitative analytical separations, and clarifying liquids. They are the standard grade cellulose filter papers for medium flow rates.

Grade 1 Qualitative Filter Papers (11 µm particle retention) are available in a variety of sizes and formats including circles, sheets and pre-pleated filters.

- Qualitative analytical separations for precipitates such as lead sulfate, calcium oxalate (hot), and calcium carbonate
- Agricultural soil analysis and seed testing procedures
- Industrial food processing techniques to extract liquid or separate solid foodstuffs from associated liquid
- Air pollution monitoring by collecting atmospheric dust from airflow and photometrically measuring the stain-intensity
- Gas detection by quantitative chromogenesis with optical reflectance
- Medium flow rate and retention for a wide range of qualitative analytical separations

**SAVE!**

Description	Code	Price	Sale
Circles, 15mm Diam [500 items]	WHAT1001-0155	<del>\$44.54</del>	<b>\$39</b>
Circles, 20mm Diam [400 items]	WHAT1001-020	<del>\$37.20</del>	<b>\$32</b>
Circles, 25mm Diam [400 items]	WHAT1001-025		
Circles, 42.5mm Diam [100 items]	WHAT1001-042		
Circles, 47mm Diam [100 items]	WHAT1001-047	<b>\$10.02</b>	<b>\$9</b>
Circles, 55mm Diam [100 items]	WHAT1001-055		
Circles, 70mm Diam [100 items]	WHAT1001-070		
Circles, 85mm Diam [100 items]	WHAT1001-085	<b>\$12.72</b>	<b>\$11</b>
Circles, 90mm Diam [100 items]	WHAT1001-090	<b>\$10.56</b>	<b>\$9</b>
Circles, 110mm Diam [100 items]	WHAT1001-110	<b>\$11.21</b>	<b>\$10</b>
Circles, 125mm Diam [100 items]	WHAT1001-125	<b>\$12.40</b>	<b>\$11</b>
Circles, 150mm Diam [100 items]	WHAT1001-150	<b>\$16.50</b>	<b>\$14</b>
Circles, 185mm Diam [100 items]	WHAT1001-185	<b>\$24.37</b>	<b>\$21</b>
Circles, 240mm Diam [100 items]	WHAT1001-240	<b>\$41.41</b>	<b>\$37</b>
Circles, 320mm Diam [100 items]	WHAT1001-320	<b>\$65.34</b>	<b>\$57</b>
Sheets, 26x31mm [1000 items]	WHAT1001-813	<b>\$77.09</b>	<b>\$67</b>
Sheets, 460x570mm [100 items]	WHAT1001-917	<b>\$136.50</b>	<b>\$118</b>



Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



**bio-strategy**  
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
31 March 2022

Ref. Code:  
WHAT-GRADE1

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.