



SOURCE

For more information and product codes please go to www.bio-strategy.com/promotions

New!

1



BioEcho Reaction Kits

Save 30% On BioEcho 50 Reaction Kits

Next Generation kits DNA/RNA extraction and analysis towards more convenience, quality and speed.

Pure DNA by only one centrifugation step.



Description	Code	Price	Sale
Blood DNA HiYield Kit (50)	BIOE010-011-050	\$260	\$182
DNA Micro Kit (50)	BIOE010-001-050	\$240	\$168
Cell Culture DNA Kit (50)	BIOE010-006-050	\$260	\$182
Plant DNA Kit (50)	BIOE010-003-050	\$300	\$210
Tissue DNA Micro Kit (50)	BIOE010-002-050	\$310	\$217
Viral RNA/DNA Swab Kit (50)	BIOE012-001-050	\$330	\$231
Viral RNA/DNA Swab Kit Plus (50)	BIOE 012-002-050	\$370	\$259

[more details](#)

[only while stocks last!](#)

In Stock!

2



PHCbi Cold Storage

Great Prices - Only While Stocks Last!



SANYMPR-N450FSH-PE
Pharmaceutical
Refrigerator / Freezer

Price \$7,795
Sale \$4,500



SANYMPR-722-PE
Pharmaceutical
Refrigerator

Price \$8,720
Sale \$5,500



SANYMDF-DU502VH-PE
Pharmaceutical
Refrigerator

Price \$18,195
Sale \$12,990



SANYMPR-S150H-PE
ECO Pharmaceutical
Refrigerator

Price \$4,360
Sale \$2,500

[View here for full specifications](#)

Clearance Sale!

3



All Stock Must Go!

For more great deals
and full list of clearance stock visit
www.bio-strategy.com/promotions



Description	Code	Price	Sale
Hei-Plate Silver 1 Package	HEID505-30080-00	\$3,477	\$1,199
Hei-Torque Core	HEID501-60410-00	\$1,667	\$929
RotaVAC Vario Control	HEID591-00141-00	\$9,396	\$4,899

[more details](#)

Best Seller!

4



Tomy Top Loading Autoclave, 79 Litre

1 Free Preventative Maintenance Visit (metro areas only)

The Tomy SX-700E top loading autoclave has an integrated (but removable) condensate bottle to trap condensate from chamber during venting.

- Mobile unit requiring no plumbing!
- Simple lid opening using either hand or foot
- Cooling fan to reduce run time
- Digital temperature & time display with analogue gauge for pressure
- Five sterilization profiles to select from (in addition to normal sterilization)
- 2-year warranty



Description	Code	Price	Sale
Autoclave, top loading, 79L	TOMYSX-700E	\$15,514	\$9,995

[more details](#)

SOURCE

For more information and product codes please go to www.bio-strategy.com/promotions

Stock Clearance!

5

grant

Grant Water Baths

Quality, Reliability & Affordability

A cost-effective range of multi-purpose water baths combining Grant's world-renowned reliability, quality and consistent performance for a wide range of laboratory applications.

- Temperature range ambient +5°C to 99°C operation
- Unique Set and Forget® technology – fast heat-up, accurate temperature control



Description	Code	Price	Sale
18L Includes clear lid, drain & base tray	GRANSAP18	\$1,809	\$1,061

[more details](#)[only while stocks last!](#)

Big Discount!

6

phcbi

PHCbi MCO-80IC CO₂ Reach-In Incubator

Specifically Designed for Critical Applications

- Copper-enriched stainless-steel germicidal interior, shelves and plenum
- Infrared CO₂ sensor
- CO₂ Range: 0 to 20%
- Field reversible door
- Accommodates roller bottle devices, 5 bottles wide x 7 bottles high

**2 In Stock**

Description	Code	Price	Sale
MCO-80IC CO ₂ Reach-In Incubator, 851 L	SANYMCO-80IC-PE	\$24,880	\$12,500

[more details](#)[price for current stock only](#)

10% Off!

7

Miele

Miele Professional PG8504 Washer

More Affordable Washing Solutions

Laboratory Underbench/Free-Standing Washer with stainless steel interior combining powerful pump with integrated heating elements.

- Door dispensers for powder detergent and neutraliser
- Suitable for operation on 240V 50Hz 13 amp supply



Description	Code	Price	Sale
PG8504 Laboratory Washer	MIEL62850422	\$7,629	\$6,865

[more details](#)

25% Off!

8


Gerhardt
Analytical Systems

Automated Kjeldahl Analysis

Expand capacity, increase throughput, and improve safety

From manual to partial-automated workflows for Kjeldahl Analysis:

- Free up time for your analysts to do more important tasks with a cost-effective system
- Improve EHS conditions with increased safety from the system design and features



Description	Price	Sale
Gerhardt Kjeldatherm KT 8s - digestion step	\$10,980	\$7,999
Gerhardt Vapodest-300 - distillation step	\$16,257	\$12,900

[more details](#)

Save!

9


Sherwood

Sherwood Scientific 926 Mk III Chloride Analyser

Chloride Analyser for the Food and Dairy Industry

The ability to read directly, salt concentration and the automatic endpoint detection system has made the 926 the instrument of choice for food manufacturers and analysts throughout the world.

- Robust and versatile coulometric method
- Fast – results in less than 40 seconds! (to 200mg/l)
- Low cost per test
- Clinical version (M926S) also available



Description	Code	Price	Sale
Model 926 Mk III Chloride Analyser	SHEW92600009	\$11,204	\$9,775

[more details](#)

Great Prices!

10

CORNING

Corning PYREX® Squat Beakers

Made from Type 1, Class A low expansion borosilicate glass, PYREX glassware has become the accepted standard in chemistry labs across the globe.

- Uniform wall thickness
- Extra large marking spot



Description	Code	Price	Sale
Beaker Squat Glass 250ml	CORN1000-250_U	\$5.80	\$3
Beaker Squat Glass 600ml	CORN1000-600_U	\$7.70	\$4
Beaker Squat H/D Glass 150ml	CORN1003-150_U	\$4.10	\$3

[more details](#)


For more information and product codes please go to www.bio-strategy.com/promotions**SOURCE****On Sale!**

11

**Adhesion Microscope Slides - In Stock**

Superfrost Polysine® electrostatically charged adhesion slides with excellent adhesive properties for high reliability in instrument-based applications.

- Adhesive coating attracts paraffin-embedded sections to slide surface
- Polysine slides are coated with a derivative of the chemistry for poly-L-lysine for tissue slices 2-5µm thick



Description	Code	Price	Sale
Slide Superfrost Polysine Coated White 72pcs	EPBRSF41296PL	\$28	\$24

[more details](#)**Great Price!**

12

**Ansell Gloves**

MICROFLEX® lightweight, yet strong nitrile glove that supports musculoskeletal health for increased worker productivity.

Only \$32!

XS-L Box 250
XL Box 230



MICRO-TOUCH® NITRA-TEX™
Excellent resistance to a broad range of commonly used chemicals. Order sizes XS - XL

Only \$19!

Box 100

Description	Code	Price	Sale
MICROFLEX Nitrile Blue XS - XL	ANSE93-833	\$47.70	\$32
Nitrile NITRA-TEX Powder-free (Medium)	ANSE4472	\$31.50	\$19

[more details](#)[price for current stock only](#)**Big Saving!**

13

**Corning Centrifuge Tubes - In Stock**

Clear polypropylene (PP) centrifuge tubes with the advanced polyethylene CentriStar™ caps.

- Superior performance
- Bulk packed 50 per bag
- Large white marking spot
- Clear accurate graduations
- Temperature range from -80°C to 120°C
- 5-year shelf life



Description	Code	Price	Sale
Tubes 15ml CentriStar [500 items]	CORN430791	\$118.80	\$84
Tubes 50ml CentriStar [500 items]	CORN430829	\$152.50	\$108

[more details](#)**In Stock!**

14

**LabAdvantage Tips**

Certified Free of:

- Pyrogens
- RNase, DNase
- PCR Inhibitors
- Human DNA



Description	Code	Price	Sale
200µl Standard Tip - Sterile Graduated Y [960 Items]	SORELAC15374	\$53.30	\$43.00
200µl Standard Tip - Graduated Natural [1000 Items]	SORELAC15724	\$19.10	\$15.90
10µl Barrier Tip - Extended Length Sterile [960 Items]	SORELAC30364T	\$107.60	\$85.30

[more details](#)**Big Discount!**

15

**DWK Laboratory Bottles**

Commonly used sizes in stock!

DWK Pyrex® Borosilicate Media Bottles, standard GL45 screw thread.

- Printed with trace code for downloadable batch COA
- Complies with ISO 4796
- For use up to 140°C



Description	Code	Price	Sale
250ml with PP Cap And Pour Ring [10 Items]	BIBB1516/06D	\$107.20	\$61
500ml with PP Cap And Pour Ring [10 Items]	BIBB1516/08D	\$120.40	\$73
1000ml with PP Cap And Pour Ring [10 Items]	BIBB1516/10D	\$169.10	\$88

[more details](#)**Buy More & Save!**

16

**DWK Glass Pasteur Pipette**

Extra savings when purchasing 5 or more packs

Disposable Kimble soda-lime glass Pasteur Pipettes for the transfer of small volumes.

- Choose from 5.75" (146mm) or 9" (228mm) plugged or unplugged
- Borosilicate pasteur pipettes also available



Description	Code	Price	Sale
146mm 5.75" Unplugged Soda-lime [1000 / pack]	KIML63A54	\$125.80	\$99
228mm 9" Unplugged Soda-lime [1000 / pack]	KIML63B93	\$154.10	\$119
228mm 9" Plugged Soda-lime [1000 / pack]	KIML63B93P	\$265.10	\$199

[more details](#)

**In Stock!**

17

**Cytiva Whatman Standard Filter Paper**

The most widely used filter paper for routine laboratory applications, qualitative analytical separations, and clarifying liquids.

They are the standard grade cellulose filter papers with low to fast flow rates.



Description		Code	Sale	Online
Grade 91	150mm Ø [1000 ITEMS]	WHAT1091-150	\$148	\$140
Grade 540	125mm Ø [100 ITEMS]	WHAT1540-125	\$45	\$43
Grade 1	240mm Ø [100 ITEMS]	WHAT1001-240	\$40	\$37
Grade 3	185mm Ø [100 ITEMS]	WHAT1003-185	\$39	\$37

[more details](#)**In Stock!**

18


METTLER TOLEDO

The Work Horse Analytical Balance**Smart Features Simplify Daily Tasks**

The touchscreen interface, easy-access menus, large digits, and results graphs simplify daily tasks. Operable with cotton, silicon, and rubber gloves.

Intuitive operation and useful features to boost efficiency.

- Accurate and reliable results
- 220 g capacity; 0.1 mg readability
- Extra large graphical color touchscreen
- Robust, full aluminum housing
- Easy cleaning



Description	Code	Sale	Online
MS204TS Analytical Balance with Touchscreen	METB30133523	\$5,665	\$5,325

[more details](#)**In Stock!**

19

**IKA Equipment**

Save up to **54%** on key items:

- Overhead Stirrers
- Homogeniser Packages
- Viscometers
- Incubator Shakers
- Popular RCT basic



Description	Code	Sale	Online
HS7 Digital Package	IKAA9026100	\$875	\$787
Stirrer Hotplate MAG RCT Basic IKAMAG®	IKAA20002620	\$850	\$765
Stirrer Magnetic C-MAG HS 7 IKAMAG®	IKAA3581200	\$710	\$639
Shaker Vortex Genius 3	IKAA3340000	\$512	\$460

[more details](#)**Can you Crack the Code?**

_____ / _____ / _____ / _____
 1 / 2 / 1 7 / 2 / 1 12 / 10 / 2 'Secret Word'
 Find it on our
 Webshop Focus
 panels

Key


First Digit Source Panel Number

Second Digit Row number (row 1 = purple header line)

Third Digit Number of word in the row

Register your answer to be in the draw to win a Food Hamper.

Entries close 17th June 2022



Register your answer at here:

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



linkedin.com/company/bio-strategy



Buy Online



View Promotions

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.