

## Corning Squat Beakers

### PYREX® Beaker with Spout

For 100 years, Corning has developed special glass for use in both chemical and life science laboratories, including PYREX® glass. Made from Type 1, Class A low expansion borosilicate glass, PYREX glassware has become the accepted standard in chemistry labs across the globe.

**Great  
Prices!**

#### Low Form, Double Scale, Graduated

Manufactured with:

- Uniform wall thickness,
- Optimum balance between thermal shock resistance and mechanical strength
- Double graduated metric scale to indicate approximate content
- Extra large marking spot

#### Heavy Duty Beakers

Specifically designed to offer:

- Best mechanical strength under harsh conditions, such as mechanized washing operations
- The 250 through 4000 ml beakers have a double graduated metric scale to indicate their approximate content.
- All sizes have an extra large marking spot

Description	Code	Price	Sale
BEAKER SQUAT GLASS 250ML [1 ITEM]	CORN1000-250_U	<b>\$5.80</b>	<b>\$3</b>
BEAKER SQUAT GLASS 600ML [1 ITEM]	CORN1000-600_U	<b>\$7.70</b>	<b>\$4</b>
BEAKER SQUAT HEAVY DUTY GLASS 150ml [1 ITEM]	CORN1003-150_U	<b>\$4.10</b>	<b>\$3</b>
BEAKER SQUAT HEAVY DUTY GLASS 600ml [1 ITEM]	CORN1003-600_U	<b>\$10.70</b>	<b>\$5</b>
BEAKER SQUAT HEAVY DUTY GLASS 1L [1 ITEM]	CORN1003-1L_U	<b>\$20.40</b>	<b>\$9</b>



Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453

E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

[shop.bio-strategy.com](http://shop.bio-strategy.com)



**bio-strategy**  
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 June 2022

Ref. Code:  
CORN-Squat

*only while stock lasts*

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.