

## Corning® Centrifuge Tube

### Superior Performance

Clear polypropylene (PP) conical bottom centrifuge tubes with the advanced polyethylene CentriStar™ caps, 25 per bulk bag.

The CentriStar™ cap has an easy-on/easy-off flat top and offers advanced ergonomics with its wider knurls and roll-over edge design for easier gripping.

- Temperature range from -80°C to 120°C
- 5-year shelf life
- Ease of Use
- Leak proof design
- Easy-on/easy-off cap
- Black printed graduations
- Large white marking spot

When used under the recommended conditions the leak-proof design, with a revolutionary plug feature, virtually eliminates all seepage.

- Reliable Quality
- RNase-/DNase-free
- Heavy-metal-free color concentrate
- Medical grade polypropylene
- Nonpyrogenic
- Sterile

**Big Savings!**

Description	Code	Price	Sale
50ml Conical CentriStar Cap [500 items]	CORN430829	<b>\$152.50</b>	<b>\$108</b>



Description	Code	Price	Sale
15ml CentriStar Cap [500 items]	CORN430791	<b>\$118.80</b>	<b>\$84</b>

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications

Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



**bio-strategy**  
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 June 2022

Ref. Code:  
CORN-Tubes

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.