

Kimble Soda-Lime Plugged and Non-Plugged Pasteur Pipettes

Disposable glass pasteur pipettes for the transfer of small volumes.

Buy More and Save!

- Choose from 5.75" (146mm) or 9" (228mm) plugged or unplugged
- Soda-Lime pipettes are manufactured from 90 expansion soda lime glass conforming to USP Type III requirements
- Borosilicate pipets are manufactured from 51 expansion, low extractable borosilicate glass conforming to USP Type 1 and ASTM E438
- 63A54P and 63B96P are cotton-plugged

Purchase 5+ Packs and Save More!



Description	Pack Size	Code	Price	Sale	5+ Price
146mm 5.75" Unplugged Soda-Lime	1000 Items	KIML63A54	\$125	\$99	\$89
146mm 5.75" Plugged Soda-Lime	1000 Items	KIML63A54P	\$239	\$199	\$179
228mm 9" Unplugged Soda-Lime	1000 Items	KIML63B93	\$154	\$119	\$109
228mm 9" Plugged Soda-Lime	1000 Items	KIML63B93P	\$265	\$199	\$179
146mm 5.75" Unplugged Borosilicate Glass Type 51	250 Items	KIML883350-0575	\$54	\$39	\$35
228mm 9" Unplugged Borosilicate Glass Type 51	250 Items	KIML883350-0009	\$59	\$39	\$35

pricing for current stock only

Contact us for more information and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2022

Ref. Code:
KIML-PIPE

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.