

Epredia Adhesion Microscope Slides


In Stock

Superfrost Polysine©

Electrostatically charged adhesion slides with excellent adhesive properties for high reliability in instrument-based applications.

- Adhesive coating attracts paraffin-embedded sections to slide surface
- Polysine slides are coated with a derivative of the chemistry for poly-l-lysine. For tissue slices 2-5µm thick

15% Off!



Description	Code	Price	Sale
Slide Superfrost Polysine Coated White 72pcs	EPBRSF41296PL	\$28	\$24

Also available Slide Printers



Thermal Transfer Microscope Slide Printer

- Efficient performance:
dual-hopper design with intelligent slide selection
- High throughput: 4-5 seconds per slide
- Simple integration:
Windows based driver and Gateway functions

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com