

Mettler Toledo MS204TS Analytical Balance

The Work Horse Analytical Balance

Smart Features Simplify Daily Tasks

The renowned MonoBloc weighing cell, FACT automatic internal adjustment, and the weighing environment filter ensure consistently reliable results.

Extra Large Graphical Color Touchscreen


The touchscreen interface, easy-access menus, large digits, and results graphs simplify daily tasks. Operable with cotton, silicon, and rubber gloves

Robust, Full Aluminum Housing

The full metal housing protects the weighing cell from impacts and environmental influences, and is resistant to harsh chemicals, including acetone.

Intuitive operation and useful features to boost efficiency

- Accurate and reliable results
- 220 g capacity; 0.1 mg readability
- Easy cleaning
- 11 applications
- EasyDirect compatible; bluetooth and wireless options



Description	Code	Price	Sale	Online
MS204TS Analytical Balance with Touchscreen	METB30133523	\$6,295	\$5,665	\$5,325

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2022

Ref. Code:
METT-Bal

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.