

LabAdvantage Tips

Standard Tip Features

- Bonded Polymer Technology
- Creates a hydrophobic surface that reduces surface tension of interior wall of the pipette tip
- Minimal sample binding
- Improves accuracy and precision
- No additives such as silicone are used
- No leaching preventing subsequent sample degradation

Barrier Tip Features:

- Small micron pore size, <20 microns
- Special membrane that prevents aerosol penetration
- Thin membrane that allows for adequate air flow to ensure accuracy and precision are not compromise

Certified Free of:

- RNase, DNase
- Human DNA
- Pyrogens
- PCR Inhibitors

**In
Stock**



Description		Code	Price	Sale
200µl Standard Tip - Sterile Graduated Yellow	[960 items]	SORELAC15374	\$53.30	\$43.00
200µl Standard Tip - Graduated Natural	[1000 items]	SORELAC15724	\$19.10	\$15.90
1000µl Standard Tip - Blue	[1000 items]	SORELAC1707	\$15.30	\$12.60
10µl Barrier Tip - Extended Length Sterile	[960 items]	SORELAC30364T	\$107.60	\$85.30

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com