

Tomy SX-700E Top Loading Autoclave

Includes two stainless steel baskets and a 2-year warranty!

The Tomy SX-700E top loading autoclave features three operating temperature ranges:

- **Warming** 45 to 95°C,
- **Heating** 45 to 104°C
- **Sterilizing** 105 to 135°C.

Safety devices include: water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection & safety valve.

- Simple lid opening using either hand or foot
- Integrated (but removable) condensate bottle to trap condensate from chamber during venting
- Cooling fan to reduce run time
- Digital temperature & time display with analogue gauge for pressure
- Five sterilization profiles to select from (in addition to normal sterilization)
- 1 **Free** Preventative Maintenance visit (*metro areas only*)

In Stock!



Description	Code	Price	Sale
Mobile 79 Litre Toploading Autoclave	TOMYSX-700E	\$15,514	\$9,995

Contact us for more information and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com