

Epredia Superfrost Adhesion Microscope Slides

Designed to secure tissue sections without the need for laboratory-applied coatings.

- Uniform coating and minimal sample loss
- For tissue slices 2-5µm thick
- Eliminates background staining in standard H&E stains



Description	Code	Price	Sale
Slide Superfrost Plus White 72pcs	EPBRSF41296SP	\$32.00	\$26.50
Slide Superfrost Polysine Coated White 72pcs	EPBRSF41296PL	\$28.00	\$24.00
Slide Superfrost Gold 72pcs	EPBRSF41296SPG	\$48.00	\$40.00

POWERED BY

SHANDON

MICROM

MENZEL-GLÄSER



Richard-Allan
Scientific

On Sale!



Enhancing precision
cancer diagnostics

The Superfrost microscope slides work best with the Epredia SlideMate Pro or SlideMate AS.



Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 September 2022

Ref. Code:
EPBR-SLIDE

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.