

Bio-Strategy Welcomes Labogene Freeze Dryers to its Portfolio!

Ask our helpful Bio-Strategy Sales Team for a configuration!

Labogene's Freeze Dryers Offer:

- Lowest condenser temperatures of -55°C, -95°C, -100°C or -110°C
- Seamless stainless-steel condenser with external cooling coils, heavily insulated to conserve energy and increase performance
- Easy draining with built-in valve, allowing removal of water after a completed freeze-drying session
- Easy operation with a built-in vacuum valve connecting the condenser and vacuum pump and allowing the condenser to cool down and the vacuum pump to warm up separately, ensuring an efficient start-up sequence

Labogene's 4-15L freeze dryers are available in 4 variants: Basic, Pro, Touch and Superior XS/XL, designed to meet the demands of laboratories of all types.

**NEW at
Bio-Strategy**



Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Rev. July 2022

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.