

## Miele Professional PG8583 ADP Washer/Disinfector

### Improving safety in your wash area.

Laboratory underbench/free-standing washer with stainless steel interior combining powerful variable speed pump with integrated heating elements enabling safe washing cycles up to 93°C.

- Integrated liquid dispensers for detergent and neutraliser
- “Touch on Steel” digital controller
- Two wash levels for greater wash load flexibility
- Underbench or standalone model
- 14 standard programs + 2 vacant program positions to meet your wash requirements
- Ecodry passive drying function – when cycle is finishing unit will automatically open door slightly at 70°C allowing residual moisture to escape easily from the wash cabinet drying the wash load
- Touch on steel controller with 3-line display for easy viewing of washer status
- Compact: 600mm x 600mm x 835mm W x D x H
- Stainless steel interior/exterior for durability and cleanliness
- Many options and accessories to select from - contact Bio-Strategy for assistance in configuring your washer

Miele offers enhanced laboratory safety by removing the manual washing step from your wash process. There is reduced likelihood of glass breakage that can result in cuts to laboratory staff. Miele is good for occupational health and safety.

**10%  
OFF!**



*Contact Bio-Strategy for optional accessories to compliment the washer.*

Description	Code	Price	Sale
Miele PG8583 ADP Undercounter/Stand Alone Laboratory Washer/Disinfector	MIEL62858306	\$10,369	\$9,332

Contact us for more information and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453  
E: sales.au@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 September 2022

Ref. Code:  
MIEL62858306

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.