

## Pall Nanosep and Microsep Centrifugal Devices

### Nanosep® Centrifugal Filters

Simple, reliable concentrating and desalting of 50 to 500µl samples.

- Rapidly concentrates 50-500µl sample volumes to 15µl
- Achieve 50X concentration and > 90% recovery in just minutes
- Available with low proteinbinding Omega™
- A wide range of MWCOs, colour-coded for easy identification
- Constructed of low-binding polypropylene
- Ultrasonically welded seals prevent bypass or seal failure
- Fits standard centrifuge rotors that accept 1.5ml tubes

Take us  
for a spin!



Description	Code	Price	Sale
Nanosep OMEGA 3k 24pk	PALLOD003C33	\$135	\$119
Nanosep OMEGA 3k 100pk	PALLOD003C34	\$453	\$400
Nanosep OMEGA 10k 24pk	PALLOD010C33	\$142	\$125
Nanosep OMEGA 10k 100pk	PALLOD010C34	\$608	\$536
Nanosep OMEGA 30k 24pk	PALLOD030C33	\$135	\$119
Nanosep OMEGA 30k 100pk	PALLOD030C34	\$453	\$400
Nanosep OMEGA 100k 24pk	PALLOD100C33	\$136	\$120
Nanosep OMEGA 100k 100pk	PALLOD100C34	\$453	\$400
Nanosep OMEGA 300k 24pk	PALLOD300C33	\$136	\$120
Nanosep OMEGA 300k 100pk	PALLOD300C34	\$453	\$400

Over  
10% Off



### Microsep™ Centrifugal Filters

Precise, quick recovery of microliter volumes.

- Rapidly concentrates 0.5-5ml sample volumes to 65–100µl
- Achieve 50X concentration and > 90% recovery in just minutes
- Deadstop feature prevents samples from spinning to dryness
- Versatile Omega™ membrane is available in a variety of MWCOs
- Colour-coded and laser etched for easy identification

Description	Code	Price	Sale
Microsep ADV OMEGA 1k 24pk	PALLMCP001C41	\$180	\$159
Microsep ADV OMEGA 1k 100pk	PALLMCP001C46	\$646	\$570
Microsep ADV Omega 3k 24pk	PALLMCP003C41	\$178	\$157
Microsep ADV Omega 3k 100pk	PALLMCP003C46	\$638	\$563
Microsep ADV Omega 10k 24pk	PALLMCP010C41	\$178	\$157
Microsep ADV Omega 10k 100pk	PALLMCP010C46	\$638	\$563
Microsep ADV Omega 30k 24pk	PALLMCP030C41	\$178	\$157
Microsep ADV Omega 30k 100pk	PALLMCP030C46	\$638	\$563
Microsep ADV Omega 100k 24pk	PALLMCP100C41	\$178	\$157
Microsep ADV Omega 100k 100pk	PALLMCP100C46	\$638	\$563

Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 September 2022

Ref. Code:  
PALL-CENT

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.