

Priorclave Q63/EH400 Front-Loading 400L Electrically Heated Autoclave

***Only 1 available at this special price!**

With cylindrical chamber 630mm Ø x 1,200mm D, the unit is suitable for a variety of loads and supplied with two shelves. Includes water cooled heat exchanger, wandering probe, printer and 5 program memory.

- Media Warming Facility
- 10-year door and chamber warranty
- Elements in chamber for ease of servicing
- Constructed to British Standards and complying with AS/NZS standards
- Grade 316L stainless steel chamber and is inherently safe with thermal and pressure interlocks
- USB port, fan forced cooling and timed free-steaming
- Single handed lever action door mechanism
- 6-8 weeks to Australia ex New Zealand warehouse

**20%
Off!***



Description	Code	Price	Sale
Autoclave, 400L Cylindrical Chamber	PRIOPS/Q63/EH400	\$49,015	
Export Crate	PRIOPC/AEC/020	\$2,217	
Wandering Probe	PRIOPC/LPT/000	\$832	
Condenser-Heat Exchanger	PRIOPC/DRC/010	\$1,933	
Tactrol 5 Program Memory	PRIOPC/PRG/TAC	\$1,652	
Printer	PRIOPC/PRN/000	\$2,004	
	Total Price:	\$57,653	\$46,122

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 September 2022

Ref. Code:
PRIO-Q63/EH400

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.