

## PHCbi LPR-400 Laboratory Refrigerator

**Only available online and while stocks last!**

Refrigerator that is designed for general laboratory storage with user-friendly features. The combination of temperature control accuracy, interior temperature uniformity, quick recovery and resistance to high ambient temperature delivers a quality storage refrigerator that characterises our commitment to engineering, storage safety and reliability.

Features:

- Cost effective model with laboratory aesthetics
- Forced air circulation technology = uniform temperature and quick temperature recovery from frequent door openings
- Condensation prevention technology
- Visual and audible alarms
- Timed cyclic defrost

Description	Code	Price	Sale
PHCbi 400L Laboratory Refrigerator	SANYLPR-400-PE	\$2,650	\$1,995

### Temperature Uniformity

Forced air circulation technology ensures a uniform temperature throughout the chamber and quick temperature recovery even with frequent door openings.

### Timer-Controlled Off-Cycle Defrosting

Defrosting starts in off mode when the accumulated compressor operation time reaches a certain threshold. During defrosting, the "dF" display and the current in-chamber temperature indications will alternate on the control panel.

### Applications

Manufactured and distributed by PHCbi, the general-purpose Laboratory Refrigerator with front glass door for easy management meets various requirement of stored goods such as cell culture media, reagents, chromatography, clinical testing kits, breast milk, medication, industrial testing materials and antibody dilutions.



**25%  
OFF!**



Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 September 2022  
or while stocks last

Ref. Code:  
SANY-LPR

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.