

Tomy Top Loading Autoclaves, 79 Litre with 2 Baskets or 36 Litre with 1 Basket

1 Free Preventative Maintenance visit (metro areas only).

The Tomy SX series top loading autoclaves are mobile and do not require plumbing. The SX-700E runs on a 230V 50Hz 15 amp circuit whilst the SX-300E runs on a 230V 50Hz 10 amp circuit.

- Three operating temperature ranges: Warming 45 to 95°C, Heating 45 to 104°C and Sterilising 105 to 135°C
- Safety devices include: water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection and safety valve
- Digital temperature and time display with analogue gauge for pressure
- Five sterilisation profiles to select from (in addition to normal sterilisation)



**40%
OFF!**



Description	Code	Price	Sale
Autoclave, Top Loading, 79 Litre	TOMYSX-700E	\$17,086	\$9,995
Autoclave, Top Loading, 36 Litre	TOMYSX-300E	\$15,929	\$9,250
Preventative Maintenance Visit	AGRMNT-S		No Charge

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 September 2022

Ref. Code:
TOMY-AUTO

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.