

NEW – Incubator Shakers

The first incubator shaker with a
patented removable shaker platform

Only \$9,995

Usually \$11,435

Code	Size
IKAA20117226	20mm
IKAA20117235	25mm

- Ø20mm or Ø25mm Shaking Stroke options
- Small Footprint
- Superior temperature performance
($\pm 0.1^{\circ}\text{C}$ stability, $\pm 0.2^{\circ}\text{C}$ uniformity)
- Temp range: Room temp+8°C to 80°C
- Decontamination (100°C*/120°C)
*When shaking platform connected




Accessories: available at additional cost



Wire grid tray



Shaker platform



Stacking kit



Roller box



Flask fixing clips



Ref. Code:
IKAA-Incu

VALID UNTIL:
30 Nov 2023



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.