

Japanese manufactured Mobile Top-loading 58L Autoclave: Quality at an affordable price!

Only \$9,995

Usually \$14,726

Code: TOMYSX-500E

The Tomy SX-500E top loading autoclave is powered from a 240V 50Hz 15amp supply and requires no plumbing so can be moved to where it is required!.


Supplied with two stainless steel baskets and a 2-year warranty!



- Digital temperature & time display with analogue gauge for pressure
- Integrated (removable) condensate bottle to trap condensate from chamber during venting
- Basket Dimensions: 300mm Ø x 181mm H
- Simple lid opening using either hand or foot
- Grade 304 stainless steel for easy cleaning and durability
- Variable speed exhaust
- Three operating temperature ranges: Warming 45 to 95°C, Heating 45 to 104°C and Sterilizing 105 to 135°C
- Safety devices include: water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection & safety valve
- Digital temperature & time display with analogue gauge for pressure
- Five sterilization profiles to select from (in addition to normal sterilization)
- Various accessories also available

Cooling Fan System

The SX series of Autoclaves are equipped as standard with a dual-fan cooling system, shortening overall operation time allowing you to focus on other things.



Ref. Code:
TOMY-SX500E

VALID UNTIL:
30 Nov 2023

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.