

Versafreeze VF15085 -86°C Underbench ULT 129L

Was \$22,760 **Now \$18,320**

Code: LAUDL003344

Laboratory upright freezer with reliable monitoring and highly efficient insulation for long-term storage and quality assurance



Features an intuitive, modern touch display with a wide range of setting options, password-protected access rights and secure monitoring via the LAUDA Cloud.

Can be built under a work bench (can easily be integrated into a laboratory line) or free-standing

Energy-saving, maintenance-free refrigeration system with fully hermetic high-performance compressors


Environmentally friendly, future-proof natural refrigerants

Protection against unauthorized access to setpoints through password-protected user management

Powerful insulation technology with high-quality materials



Working temperature
min. -86 °C
max. -50 °C



Insulation (right to left)
Interior with stainless steel walls, tightly meshed evaporator, thermal film, vacuum panels, impermeable foam, appliance front



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
LAUD-Versa

VALID UNTIL:
15 Dec 2023