

Extreme Environment Shaker

Extreme Performance, Under Extreme Conditions!

Was **\$8,300** Now **\$5,995**


Code: OHAU30391817



watch >

The SHEX1619DG orbital shaker's unique construction reliably delivers exceptional speed control, accuracy and durability under the most extreme environmental conditions.

- Motion: Orbital, 19 mm
- Capacity: 16 kg
- Speed Range: 15 rpm – 500 rpm
- Tray size: 33 cm x 33 cm
- Timer: 1 second – 160 hours
- Includes non-slip rubber mat



The extended ribbon cable (1.68m) can easily pass underneath an incubator door or via incubator's utility port, for external control. Avoid unnecessary disruption of the internal environment.

For use in cold rooms and incubators, from -10°C to 60°C, CO₂ environments, and up to 100% humidity

Various clamps, racks and accessories available – please contact us for more information.



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
OHAU-SHEX

VALID UNTIL:
15 Dec 2023