

## VIP ULT -86°C Chest Freezer


**Small size, small price**

Was \$9,040

**Now \$5,876**

Code: SANYMDF-C8V1-PE

*Only while stocks last*



Ideal for storage of small volume clinical trial samples and temperature control testing.

- Temp Control Range -60°C to -80°C
- Filter-less design
- Capacity 84 L (42 x 2' cryoboxes)
- 5-year warranty



The new generation vacuum panels for increased energy efficiency and cooling performance

### Related Products

Cryo boxes in a range of sizes available from various suppliers



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
SANY-84LChest

VALID UNTIL:  
15 Dec 2023