

VIP Eco ULT Freezer 528L

Preservation sub-zero



Was \$17,385 **Now \$11,300**

Code: SANYMDF-DU502VH-PE

Only while stocks last

Sustainability Savings!

Save now and into the future with the ULT -80°C Freezer

- Save on energy costs (6.7 kWh/day)
- Save on HVAC costs (419 W heat dissipation)
- Save now with our special price!

Features:

- Temp Control Range -40°C to -86°C
- Internal temperature/door/error logging touch pad interface
- Capacity 528 L (384 x 2' cryoboxes)
- 5-year warranty



Related Products

Cryo boxes and racks in a range of sizes
available from various suppliers



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
SANY-VIP_528L

VALID UNTIL:
15 Dec 2023