

Hub RN Needles

Recently purchased a Hamilton syringe with removable needle and looking for spares or want to re-stock?

Check out our special pricing across all standard small and large hub RN needles

From \$74

Removable needles give you the best of both worlds

1. the benefit of a replaceable needle for the occasional bend or clog
2. the minimal dead volume of a cemented needle. The unique removable needle design installs the needle precisely at the zero line

Large



Compatible with syringes ranging from 250uL up to 10mL.

Small



Compatible with syringes ranging from 2.5uL up to 100uL.

Need custom needles, just give us a call

see pricing list page 2

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code: **VALID UNTIL:**
HAMC-HubRN **30 Sep 2024**

Hub RN Needles

All come in a pack of 6

Code	Size	Price	Special
Large			
HAMC7779-02	26s/2"/2	105	79
HAMC7779-04	26/2"/2		
HAMC7779-03	22s/2"/2	142	106
HAMC7779-01	22/2"/2	105	79
HAMC7780-01	26s/2"/3	98	74
HAMC7780-02	26/2"/3		
HAMC7780-03	22s/2"/3	148	111
HAMC7780-04	22/2"/3	108	81
HAMC7784-03	26s/2"/5	225	169
HAMC7784-04	26/2"/5		
HAMC7784-01	22s/2"/5	262	196
HAMC7784-02	22/2"/5	225	169
Small			
HAMC7758-02	26s/2"/2	116	87
HAMC7758-04	26/2"/2		
HAMC7758-03	22s/2"/2	156	117
HAMC7758-01	22/2"/2	116	87
HAMC7762-06	33/1.5"/3	167	125
HAMC7762-05	32/2"/3		
HAMC7762-04	31/2"/3		
HAMC7762-03	30/2"/3	135	101
HAMC7762-02	28/2"/3		
HAMC7762-01	7/2"/3		
HAMC7768-01	26s/2"/3	98	74
HAMC7768-02	26/2"/3		
HAMC7770-01	22s/2"/3	135	101
HAMC7770-02	22/2"/3	98	74
HAMC7784-07	26s/2"/5	225	169
HAMC7784-08	26/2"/5		
HAMC7784-05	22s/2"/5	262	196
HAMC7784-06	22/2"/5	225	169



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
HAMC-HubRN

VALID UNTIL:
30 Sep 2024