

## PG8504 Laboratory Underbench/ Free-Standing Washer

**Only \$7,642**

Was \$8,491

Code: MIEL62850422  
(without accessories)

**Are you still manually  
washing glassware?**

**Let's put some precision and  
safety into this procedure!**

Miele offers enhanced laboratory safety by removing the manual washing step from your wash process. There is reduced likelihood of glass breakage that can result in cuts to lab staff.

### **Miele - good for occupational health and safety**

With stainless steel interior combining powerful 500L/ min pump with integrated heating elements enabling safe washing cycles with temperatures up to 70°C .

- Enhanced safety by removing the manual washing step from your wash process
- Door dispensers for powder detergent and neutraliser
- Suitable for operation on 240V 50Hz 13 amp supply
- 6 standard programs + 2 vacant program positions to meet your wash requirements
- Touch on steel controller with 3-line display for easy viewing of washer status



Many options and accessories to select from, contact Bio-Strategy for assistance in configuring your washer.

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
MIEL-PG8504

VALID UNTIL:  
30 Apr 2024