



NEW

FrostLess 725L -86°C Freezer

An optimal choice for ultra-low temperature storage, offering less frost build-up than our category-leading models, while maintaining the reliability and energy efficiency today's laboratories require.

Combat frost build-up damage to freezer components with this low frost model.

Frost Resistant Inner Doors

Frost accumulation may cause damage and deterioration of freezer components.

The **all new** thermal insulation design prevents large amounts of frost buildup on the outside of the inner doors reducing maintenance requirements

Less frost buildup puts less stress on the outer door and latching mechanism to help preserve optimal performance.

Frost Mitigation

Using the residual heat from the refrigeration system, the perimeter heater under the gasket allows for gentle warmth reducing ice buildup on the door gasket.

Less frost build-up reduces manual labor and improves lab workflow.



Interested?
Give us a call

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref:
Feb 2024