

Frontier™ 2000 Multi Centrifuge

**Economical, Non-refrigerated
Multipurpose Benchtop Centrifuge Package**

Only \$1,950

Was \$2,570

Code: OHAU30602527

Includes:

FC2706 Centrifuge & 12x15ml Fixed Angle Rotor



A reliable low-speed general purpose centrifuge, the FC2706 Multi Centrifuge is economical, intuitive, and designed with convenient operation in mind.

The included 15ml rotor can accommodate both round and conical bottom tubes (including common blood tubes).

- Reaches a maximum speed of 5,000 rpm/3,080 x g in 12 sec
- Unique anit-slam lid design, remains open without accidentally closing
- Automatic lid opening once the run is complete
- Backlit LCD display, allows the user to easily view settings as well as adjust speed/time during operation
- Housing, chamber and rotor buckets are constructed of corrosion-resistant stainless steel that holds up well to organic solvents.
- Imbalance sensor, for both user and instrument safety
- Versatile, with 5 x different rotors available
- 2-year warranty

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
OHAU-FC2706

VALID UNTIL:
30 Apr 2024