

## 79L Top Loading Autoclave

**Only \$10,975**

Was \$17,342

Code: TOMYSX-700E

**Supplied with 2 stainless steel baskets and a 2-year warranty!**



The Tomy SX-700E top loading autoclave is powered from a 240V 50Hz 15amp supply and requires no plumbing so can be moved to where it is required!.

- Integrated (removable) condensate bottle to trap condensate from chamber during venting
- Basket Dimensions: 345mm Ø x 181mm H
- Grade 304 stainless steel for easy cleaning and durability
- Variable speed exhaust
- Three operating temperature ranges: Warming 45 to 95°C, Heating 45 to 104°C and Sterilizing 105 to 135°C
- Safety devices include: water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection & safety valve.
- Digital temperature & time display with analogue gauge for pressure.
- Five sterilization profiles to select from (in addition to normal sterilization)
- Cooling fan to reduce run time

### Various accessories also available:

Including punched metal stainless steel baskets, solid stainless steel baskets, hybrid basket: solid lower third with punched metal stainless steel upper and waste containers



**Simple lid opening using either hand or foot**

### HybridBaskets



### Waste Containers



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
TOMY-SX-700E

VALID UNTIL:  
31 May 2024