

# Miele

## PG 8583 ADP Washer/Disinfector

**Only \$10,647**

Usually \$11,830

Code: MIEL62858306

Includes **FREE**  
Detergent & Neutraliser



- Compact: 600mm (W) x 600mm (D) x 850mm (H) (820mm with lid removed)
- Two wash levels
- Spray arm and spray pressure monitoring ensuring a clean wash
- Stainless steel interior/exterior for durability and cleanliness
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- Includes pump for AD water and waterproofing kit
- Dispensing pumps for liquid detergent and liquid neutraliser ensure accurate dosing
- Many options and accessories to select from
- 2-year Warranty

**Passive EcoDry drying function** – activates as temperature descends to 70°C opening the door slightly enabling the load to be dried

**15%** discount on optional accessories

Give us a call to discuss your complete washing solution



**bio-strategy**

Part of DKSH Group

Ref. Code:  
MIEL:PG8583

VALID UNTIL:  
31 Dec 2024



1800 00 84 53

sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.