

165L
~~Only \$6,204~~ **Now \$4,407**
Usually \$7,755

Code: SANYMCO-170AICL-PE
Ambient Temp +5°C to 50°C



230L
~~Only \$7,328~~ **Now \$5,105**
Usually \$9,160

Code: SANYMCO-230AICL-PE
Ambient Temp +5°C to 50°C ± 0.25°C



Other models
also available
Ask us about
special pricing

Hit a PAR with your cell culture!

- **Passive Contamination Control**
Copper-Embedded Stainless-Steel interior & shelves
- **Active Contamination Control**
Airflow Sterilised by UV light (Optional)
- **Reactive Contamination Control**
Vapourised H₂O₂ decontamination (Optional)

Benefits:

- Removable water pan
- Clear and large display LCD Touch Screen
- Dual IR CO₂ sensor for fastest CO₂ recovery
- Double stackable

Code	Optional Upgrades	Was	Was	Now
SANYMCO-170UVS-PE	UV System Kit	\$810	\$648	\$495
SANYMCO-170HB-PE	H ₂ O ₂ Decontamination Control Board	\$765	\$612	\$420
SANYMCO-170EL-PW	Electric Lock - H ₂ O ₂ Upgrade	\$775	\$620	\$470