

V-200 Vortex Mixer

FREE
Accessory Cup Pack
(while stock lasts)

Only \$575

Usually \$821

Code: COLE-04729-03



Compact unit perfect for individual work stations.

- Mix a wide variety of tube sizes and shapes with the included rubber cup head
- Heavy die-cast body won't dance across your bench
- Automatic start when the soft rubber cup is pressed down
- Variable speed mixing from 200 to 2500 rpm
- Select between "touch" mode or continuous operation
- 3-year warranty

FREE Accessory Cup Pack Includes:

- 1 x plastic cradle for standard microtitre plate
- 1 x circular foam insert for beakers and flasks (up to 500ml capacity)
- 1 x rectangular foam insert (holds 8 x 0.2ml, 8 x 0.5ml and 16 x 1.5ml microcentrifuge tubes)
- 1 x solid rectangular foam insert (blank for custom holes)

Please note: The microtitre plate holder SA8/1 should not be used at speeds above 800rpm

Ref. Code: COLE-V200

Valid Until: 28 Feb 2025

