



# Next Generation Moisture Analysers

With carbon fibre heating technology for fast,  
reliable, and accurate results

Moisture analysers are an excellent alternative to traditional oven drying methods and offer a user-friendly option compared to wet chemistry analysers. Whether you need a model for repetitive production testing or for method development in the lab, you can choose the one that best meets your needs!

**Only \$1,995**

**Usually \$2,800**

Code: OHAU30971934

Model: MB32



Perfect for  
routine moisture testing

- Temperature Range: 40 - 180 °C
- Drying Programs: Standard, Fast
- Shut-off Criteria: A30, A60, A90, Manual, Timed
- Method Library: 2
- Result Storage: 1
- Knob control
- 4.3" LCD display

**Only \$2,950**

**Usually \$4,100**

Code: OHAU30971935

Model: MB62



Designed for durability and  
flexible performance in a variety  
of moisture applications

- Temperature Range: 40 - 200 °C
- Drying Programs: Standard, Fast, Ramp, Step
- Shut-off Criteria: A10, A30, A60, A90, A150, Manual, Timed
- Method Library: 20
- Result Storage: 100
- Knob control
- 5" LCD display

**Only \$3,350**

**Usually \$4,650**

Code: OHAU30971936

Model: MB92



High performance moisture  
analysis with advanced features  
for the most demanding  
application

- Temperature Range: 40 - 200 °C
- Drying Programs: Standard, Fast, Ramp, Step
- Shut-off Criteria: A10, A30, A60, A90, A150, Manual, Timed, Free mgs
- Method Library: 20
- Result Storage: 200
- Colour easy-to-use Touchscreen display

Ref. Code: OHAU-MB

Valid Until: 31 Mar 2025



Part of DKSH Group



1800 00 84 53  
sales.au@bio-strategy.dksh.com  
www.bio-strategy.com  
shop.bio-strategy.com

