

### 58L

**Only \$8,975**

**Usually \$14,726**

Code: TOMYSX-500E

- Basket: 300mm Ø x 181mm H
- Chamber: Ø x H: 325 x 733mm
- External: 410 x 574 x 970mm W x D x H

### 79L

**Only \$10,795**

**Usually \$17,342**

Code: TOMYSX-700E

- Basket: 345mm Ø x 181mm H
- Chamber: Ø x H: 370 x 774mm
- External: 470 x 625 x 1,003mm W x D x H

The SX-500E & SX-700E top loading autoclaves have castors for mobility and a removable condensate bottle negating the need for drains and only require a 240V 50Hz 15amp supply.

**Both models include  
two baskets & 2-year warranty**

### Features common to both Autoclaves:

- Simple lid opening using either hand or foot
- Safety devices include: water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection & safety valve
- Five sterilization profiles to select from (in addition to normal sterilization)
- Three operating temperature ranges: Warming 45 to 95°C, Heating 45 to 104°C & Sterilizing 105 to 135°C
- Variable speed exhaust and cooling fan to reduce run time
- Grade 304 stainless steel for easy cleaning and durability
- Digital temperature & time display with analogue gauge for pressure



TOMYSX-500E



TOMYSX-700E

Ref. Code: TOMY-SX

Valid Until: 31 Mar 2025