



AQUASEARCHER™ AB33M1

Easy-to-Use & Accurate
Multi-Parameter Benchtop Meter



Only \$1,395

Usually \$1,827

Code: OHAU30589825

Includes ST310 pH electrode,
and STCON3 conductivity probe

- Two Independent Channels
- Calibration due alarm
- 1,000 measurement memory
- RS232 & USB ports as standard
- pH Measurement Range:
-2.00 – 20.00 pH
- Conductivity measuring range:
0.01 μ S/cm – 500 mS/cm

The multi-parameter bench meter designed with reliable measurement and user-interface in mind.

The i-Steward electrode-monitoring software ensures repeatable and consistent measurement accuracy for optimal peace of mind.

The 6.5 inch LCD display and touch keypad makes changing the parameters, executing setup, and performing calibration as simple as using a smartphone.

Ref. Code: OHAU-AB33M1

Valid Until: 30 Apr 2025



1800 00 84 53
sales.au@bio-strategy.dksh.com
www.bio-strategy.com
shop.bio-strategy.com



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.