

Stock Clearance

Only \$24

Usually \$30

Pack of 5

Code: AZLOBWP1000P

**Ideal for a wide range of liquid storage,
media preparation and sampling applications**



- Translucent, rigid polypropylene round bottles with polypropylene screw caps
- Leak-proof linerless cap design
- Smooth inner contours of neck for efficient pouring and ease of cleaning
- Heavy duty and durable in use
- Suitable for food contact applications in accordance with EU and U.S. FDA regulations
- Good chemical resistance
- Autoclavable at 121oC

NOTE: When autoclaving bottles always ensure the caps are completely loosened or removed to prevent accidental implosion

We offer a wide range of storage bottles to meet all your laboratory needs.
Chat to your local rep to find the right storage bottles for your requirements.

Ref. Code: DWK_WideNeck

Valid Until: 31 July 2025