

Miele

PG8583 CD Washer/Disinfector

With DryPlus HEPA filtered drying
No need for a separate oven to dry your loads!



Only \$19,995

Usually \$22,300

Code: MIEL62858351

Complete system requires optional
baskets and inserts



FREE

Detergent & Neutraliser



Plus 15%

discount on optional accessories

**Please call us to discuss your
complete washing solution**

- Conductivity monitoring to ensure load is clean
- Dispensing pumps for automatic dispensing of liquid detergent and liquid neutraliser
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- 14 standard programs + 2 freely configurable programs
- Two wash levels configured for washing wide neck and narrow neck glassware

- Spray arm & spray pressure monitoring ensuring a clean wash
- Underbench or standalone model
- Stainless steel interior/exterior for durability and cleanliness
- ADP model including pump for AD water and waterproofing kit
- Many options and accessories to select from
- Pull out drawer for housing detergent and neutraliser

Ref. Code: MIEL-PG8583CD

Valid Until: 30 June 2025

 **bio-strategy**
Part of DKSH Group



1800 00 84 53
sales.au@bio-strategy.dksh.com
www.bio-strategy.com
shop.bio-strategy.com

 **DKSH**

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.